

CCHFV

- Detection of Crimean Congo Haemorrhagic Fever Virus (CCHFV) specific RNA
- Parallel testing within tropical RT-PCR test kit panel
- Ready-to-use kit including Internal and Positive Control
- Evaluated at different reference laboratories
- Compatible with various Real-Time PCR platforms

RealStar® CCHFV RT-PCR Kit 1.0

Prod. No 181003

For use with:

ABI Prism® 7500 SDS and 7500 Fast SDS, m2000rt, Mx 3005P™ QPCR System, VERSANT™ kPCR Molecular System AD, LightCycler® 480 Instrument II, Rotor-Gene™ 3000/6000, Rotor-Gene™ Q 5/6 plex Platform

RealStar® CCHFV RT-PCR Kit 1.2



Prod. No 181202

For use with:

SmartCycler® II, LightCycler® 1.2 / 1.5 / 2.0

always a drop ahead.

RealStar® CCHFV RT-PCR Kit



Dilution series of CCHFV specific RNA, ranging from 10^1 to 10^4 copies/ μ l, detected in the FAM detection channel [A]. The Internal Control (IC) is detected in the JOE detection channel [B].

