

Clostridium difficile

- Detection and differentiation of *C. difficile* Toxin A and Toxin B specific DNA
- Parallel testing within enteric virus and bacteria test kit panel
- Ready-to-use kit including Internal and Positive Control
- Compatible with various Real-Time PCR platforms
- CE-IVD marked *in-vitro* diagnostic test




RealStar® *Clostridium difficile* PCR Kit 1.0 Prod. No 171013

For use with:

ABI Prism® 7500 SDS and 7500 Fast SDS, m2000rt, Mx 3005P™ QPCR System, VERSANT™ kPCR Molecular System AD, LightCycler® 480 Instrument II, Rotor-Gene™ 3000/6000, Rotor-Gene™ Q 5/6 plex Platform

always a drop ahead.

RealStar® *Clostridium difficile* PCR Kit



C. difficile Toxin B specific DNA detected in the FAM detection channel [A] and *C. difficile* Toxin A specific DNA detected in the Cy5 detection channel [B]. The Internal Control (IC) is detected in the JOE detection channel [C] and impaired in case of high positive samples.

Red: *C. difficile* A⁺B⁺; Blue: *C. difficile* A⁻B⁺; Green: *C. difficile* A⁻B⁻; Black: NTC.

