

PIV

- Detection and differentiation of human Parainfluenza Virus 1–4 specific RNA
- Parallel testing within respiratory RT-PCR test kit panel
- Ready-to-use kit including Internal and Positive Controls
- Compatible with various Real-Time PCR platforms
- CE-IVD marked *in-vitro* diagnostic test


RealStar® PIV RT-PCR Kit 1.0 Prod. No 261013

For use with:

ABI Prism® 7500 SDS and 7500 Fast SDS, m2000rt, Mx 3005P™ QPCR System, VERSANT™ kPCR Molecular System AD, LightCycler® 480 Instrument II, Rotor-Gene™ 3000/6000, Rotor-Gene™ Q 5/6 plex Platform

always a drop ahead.

RealStar® PIV RT-PCR Kit



Dilution series (10^0 to 10^4 copies/ μ l) of PIV Genus *Respirovirus* specific RNA shown as PIV-1 [A] and PIV-3 [B] in the FAM detection channel and PIV Genus *Rubulavirus* specific RNA, shown as PIV-2 [C] and PIV-4b [D] in the Cy5 detection channel. The Internal Control is detected in the JOE detection channel [E].

