



**WNV**



**RealStar® WNV RT-PCR Kit**

- Detection of West Nile Virus (WNV) specific RNA
- Detection of WNV Lineage 1 and WNV Lineage 2
- Ready-to-use kit including Internal and Positive Control
- Compatible with various Real-Time PCR platforms



**RealStar® WNV RT-PCR Kit 1.0**  
Prod. No 321003

**For use with:**

ABI Prism® 7500 SDS, m2000rt,  
Mx 3005P™ QPCR System,  
VERSANT™ kPCR Molecular System AD,  
LightCycler® 480 Instrument II,  
Rotor-Gene™ 3000/6000,  
Rotor-Gene™ Q 5/6 plex Platform

WNV specific RNA detected in the FAM detection channel [A]. Green: Positive Control; Red: WNV Lineage 1; Brown: WNV Lineage 2. The Internal Control (IC) is detected in the JOE detection channel [B].

