

# AltoStar® – Product Overview CE-IVD



## Sample matrices with projected launch dates\*

The table below provides an overview about the human specimen types, which are validated (✓) or when they will be validated for use with the respective AltoStar® (RT-)PCR Kit.

AltoStar® (RT-)PCR kit	Order No.	Plasma	Whole Blood	Serum	Urine	Swabs	BAL	Stool	CSF
Adenovirus 1.5	AS0301543	✓	✓						
<i>alpha</i> Herpesvirus 1.5	AS0081543					✓			
BKV 1.5	AS0031543	✓	✓		✓				Q3/Q4 2023
CMV 1.5	AS0021543	✓	✓		✓				Q3/Q4 2023
EBV 1.5	AS0131543	✓	✓						Q3/Q4 2023
HAV 1.5	AS0241543	✓							
HBV 1.5	AS0201513	✓							
HCV 1.5	AS0211513	✓		✓					
HEV 1.5	AS0271543	✓							
HHV-6 1.5	AS0311543	✓	✓						Q2/Q3 2023
HIV 1.5	AS0221513	✓							
HSV 1.5	AS0061543	✓							Q2/Q3 2023
Influenza S&T 1.5	AS0161543					✓			
JCV 1.5	AS0041543	✓							Q2/Q3 2023
MERS-CoV 1.5	AS0391543					✓			
Norovirus 1.5	AS0051543							✓	
Parvovirus B19 1.5	AS0101543	✓							
SARS-CoV-2 1.5	AS0821543, AS0821547					✓			
VZV 1.5	AS0071543	✓							Q2/Q3 2023
Pneumocystis jirovecii 1.5**	AS0551543						Q2/Q3 2023		

AltoStar®

altona Diagnostics GmbH  
Mörkenstr. 12  
22767 Hamburg  
Germany

phone +49 40 548 0676 0  
fax +49 40 548 0676 10  
e-mail [info@altona-diagnostics.com](mailto:info@altona-diagnostics.com)  
[www.altona-diagnostics.com](http://www.altona-diagnostics.com)

\*Launch dates are prone to changes!

\*\*Product not CE-IVD marked yet. CE-IVD marking planned to be achieved with the launch date specified for the respective specimen type

altona®, AltoStar® and all associated marks are trademarks of altona Diagnostics GmbH.

All other registered names and trademarks used in this document are not to be considered unprotected by law.

© altona Diagnostics GmbH 2022. All rights reserved.

LIS-AS-Matrices-EN-02\_08/2022