

EHEC

- Detection and differentiation of Shiga Toxin 1 and Shiga Toxin 2 and ipaH specific DNA
- Parallel testing within enteric virus and bacteria test kit panel
- Ready-to-use kit including Internal and Positive Control
- Compatible with various Real-Time PCR platforms
- CE-IVD marked *in-vitro* diagnostic test


RealStar® EHEC PCR Kit 2.0

Prod. No 292013

For use with:

ABI Prism® 7500 SDS and 7500 Fast SDS, Mx 3005P™ QPCR System, VERSANT™ kPCR Molecular System AD, LightCycler® 480 Instrument II, Rotor-Gene™ 3000/6000, Rotor-Gene™ Q 5/6 plex Platform, CFX96™/Dx Real-Time System

RealStar® EHEC PCR Kit



Dilution series (10^1 to 10^6 copies/ μ l) of Shiga Toxin 1, Shiga Toxin 2 and ipaH specific DNA detected in [A] the Cy5 detection channel, [B] the FAM detection channel and [C] the ROX detection channel, respectively. The Internal Control (IC) is detected in the JOE detection channel [D]. Black: NTC.

