

JCV

- Detection and quantification of JC Virus (JCV) specific DNA
- Parallel testing within Herpes- and Polyomavirus test kit panel
- Ready-to-use kit including Internal Control and Quantification Standards
- Compatible with various Real-Time PCR platforms
- CE-IVD marked *in-vitro* diagnostic test

RealStar® JCV PCR Kit 1.0

Prod. No 041013

For use with:

ABI Prism® 7500 SDS and 7500 Fast SDS,
m2000rt, Mx 3005P™ QPCR System,
VERSANT™ kPCR Molecular System AD,
LightCycler® 480 Instrument II,
Rotor-Gene™ 3000/6000,
Rotor-Gene™ Q 5/6 plex Platform

RealStar® JCV PCR Kit 1.2

Prod. No 041212



For use with:

SmartCycler® II, LightCycler® 1.2 / 1.5 / 2.0

always a drop ahead.



RealStar® JCV PCR Kit



Dilution series of JCV specific DNA, ranging from 10^0 to 10^9 copies/ μ l, detected in the FAM detection channel [A]. The Internal Control is detected in the JOE detection channel [B].

