

MERS-CoV

- Detection of Middle East Respiratory Syndrome coronavirus specific RNA
- Ready-to-use kit including Internal and Positive Control
- Compatible with various Real-Time PCR platforms
- CE-IVD marked *in-vitro* diagnostic test



RealStar® MERS-CoV RT-PCR Kit 1.0 Prod. No 391012

For use with:

ABI Prism® 7500 SDS and 7500 Fast SDS,
Mx 3005P™ QPCR System,
VERSANT™ kPCR Molecular System AD,
LightCycler® 480 Instrument II,
Rotor-Gene™ 3000/6000,
Rotor-Gene™ Q 5/6 plex Platform

always a drop ahead.

RealStar® MERS-CoV RT-PCR Kit



Amplification signals of the two different MERS-CoV targets *upE* (blue) and *Orf1a* (red) in the FAM-channel [A]. The internal control signal is detected in the JOE-channel [B].

