

Norovirus




- Detection of human pathogenic Norovirus specific RNA
- Differentiation of Norovirus Genogroup I and Norovirus Genogroup II
- Ready-to-use kit including Internal and Positive Controls
- Compatible with various Real-Time PCR platforms
- CE-IVD marked *in-vitro* diagnostic test

RealStar[®] Norovirus RT-PCR Kit 2.0 Prod. No 052013

For use with:

ABI Prism[®] 7500 SDS and 7500 Fast SDS, m2000rt, Mx 3005P[™] QPCR System, VERSANT[™] kPCR Molecular System AD, LightCycler[®] 480 Instrument II, Rotor-Gene[™] 3000/6000, Rotor-Gene[™] Q 5/6 plex Platform

RealStar[®] Norovirus RT-PCR Kit



Dilution series (10^0 to 10^4 copies/ μ l) of Norovirus genogroup II specific RNA detected in the FAM detection channel [A] and Norovirus genogroup I specific RNA detected in the Cy5 detection channel [B]. The Internal Control (IC) is detected in the JOE detection channel [C].

