

## Pneumocystis jirovecii

- Detection and quantification of *Pneumocystis jirovecii* specific DNA
- Ready-to-use kit including Internal Control and Quantification Standards
- Compatible with various Real-Time PCR platforms
- CE-IVD marked *in-vitro* diagnostic test



### RealStar® Pneumocystis jirovecii PCR Kit 1.0 Prod. No 551013

**For use with:**

ABI Prism® 7500 SDS and 7500 Fast SDS,  
Mx 3005P™ QPCR System,  
VERSANT® kPCR Molecular System AD,  
LightCycler® 480 Instrument II,  
Rotor-Gene® 6000,  
Rotor-Gene® Q 5/6 plex Platform,  
CFX96™ Real-Time System,  
CFX96™ Deep Well Real-Time System

always a drop ahead.

## RealStar® Pneumocystis jirovecii PCR Kit



Dilution series of *Pneumocystis jirovecii* specific DNA, ranging from  $10^1$  to  $10^8$  copies/ $\mu$ l, detected in the FAM detection channel [A].

The Internal Control is detected in the JOE detection channel [B].

