



**VZV**

- Detection and quantification of Varicella-Zoster Virus (VZV) specific DNA
- Parallel testing within Herpes- and Polyomavirus test kit panel
- Ready-to-use kit including Internal Control and Quantification Standards
- Compatible with various Real-Time PCR platforms
- CE-IVD marked *in-vitro* diagnostic test



**RealStar® VZV PCR Kit**



**RealStar® VZV PCR Kit 1.0**  
Prod. No 071013

**For use with:**

ABI Prism® 7500 SDS and 7500 Fast SDS, m2000rt, Mx 3005P™ QPCR System, VERSANT™ kPCR Molecular System AD, LightCycler® 480 Instrument II, Rotor-Gene™ 3000/6000, Rotor-Gene™ Q 5/6 plex Platform

Dilution series of VZV specific DNA, ranging from  $10^0$  to  $10^8$  copies/ $\mu$ l, detected in the FAM detection channel [A]. The Internal Control is detected in the JOE detection channel [B].

