



Yellow Fever Virus

- Detection of Yellow Fever Virus specific RNA
- Ready-to-use kit including Internal and Positive Control
- Compatible with various Real-Time PCR platforms
- CE-IVD marked *in-vitro* Diagnostic Test

RealStar[®] Yellow Fever Virus RT-PCR Kit



YFV specific RNA detected in the FAM detection channel [A], black: Positive Control, blue: Negative Control. The Internal Control (IC) is detected in the JOE Channel [B].

RealStar[®] Yellow Fever Virus RT-PCR Kit 1.0 Prod. No 671013

For use with:

Mx 3005P™ QPCR System,
VERSANT[®] kPCR Molecular System AD,
ABI Prism[®] 7500 SDS and 7500 Fast SDS,
LightCycler[®] 480 Instrument II,
Rotor-Gene[®] 6000,
Rotor-Gene[®] Q 5/6 plex Platform,
CFX96™ Real-Time PCR Detection System,
CFX96™ Deep Well Real-Time PCR
Detection System

always a drop ahead.

altona Diagnostics GmbH
Mörkenstr. 12
22767 Hamburg
Germany

phone +49 40 548 06 76 - 0
fax +49 40 548 06 76 - 10
e-mail info@altona-diagnostics.com
www.altona-diagnostics.com

