





alpha Herpesvirus

- Detection of human *alpha* Herpesvirus specific DNA
- Differentiation of HSV-1, HSV-2 and VZV in one tube
- Ready-to-use kit including Internal and Positive Controls
- Compatible with various Real-Time PCR platforms
- CE-IVD marked *in-vitro* diagnostic test



RealStar[®] alpha Herpesvirus PCR Kit



RealStar[®] alpha Herpesvirus PCR Kit 1.0 Prod. No 081013

For use with:

ABI Prism[®] 7500 SDS and 7500 Fast SDS, m2000rt, Mx 3005P™ QPCR System, VERSANT™ kPCR Molecular System AD, Rotor-Gene™ 3000/6000, Rotor-Gene™ Q 5/6 plex Platform, CFX96™/Dx Real-Time System

Dilution series (10^0 to 10^4 copies/ μ l) of VZV, HSV-1 and HSV-2 specific DNA detected in [A] the FAM detection channel, [B] the ROX detection channel and [C] the Cy5 detection channel, respectively. The Internal Control is detected in the JOE detection channel [D].

