



AltoStar® Automated Whole Blood Pretreatment

A software that provides a solution for the AltoStar® Workflow to perform the pretreatment of whole blood samples in an automated way.

- No need for manual pipetting
- Pipetting done automatically by the AltoStar® Automation System AM16
- Mixing of clot free whole blood and AltoStar® Whole Blood Pretreatment Buffer 1.5 by pipetting up and down
- Subsequent nucleic acid extraction possible

Reagents and consumables

AltoStar® Whole Blood Pretreatment Buffer 1.5*

Prod. no. WBPB15-06
576 reactions
Buffer for the pretreatment of whole blood samples

AltoStar® Purification Kit 1.5*

Prod. no. PK15-06
1152 reactions
Reagents for the extraction of nucleic acids

Secondary Sample Tubes with Barcode

Prod. no. VK000043
1000 tubes and lids
Prelabeled tubes with a barcode



**For further information
please contact your local
sales representative.**

altona Diagnostics GmbH
Mörkenstr. 12
22767 Hamburg
Germany

phone +49 40 548 0676 0
fax +49 40 548 0676 10
e-mail info@altona-diagnostics.com
www.altona-diagnostics.com

*Products for research use only. Not for use in diagnostic procedures.

AltoStar® and all associated marks are registered trademarks of altona Diagnostics GmbH.

© altona Diagnostics GmbH 2022. All rights reserved.

FLY-01-AS-WBPB15-06-EN-01_10/2022